

Swissedent Ultra Sonic Cleaner

Swissedent Ultra Sonic Cleaner works on the principle of sonic wave form which greatly enhances the cleaning of all types of dentures, orthodontic appliances (braces) removable bridges, gum shields, mouth guards and small jewellery.

It is vital that dentures are cleaned regularly and thoroughly to minimize the risk of infection and bad breath odours. Many dentures wearers find it is embarrassing to soak their dentures overnight. Sonic cleaning means an end to overnight soaking. Just 5 minutes twice daily using your Swissedent Sonic Denture Cleaner with the recommended denture cleaning tablets/sachet will keep your dentures or other appliances sparkling clean and fresh.

Instructions:

1. Remove the perforated tray from your Swissedent Sonic Cleaner and place your dentures inside (perforated tray is used for small items).
2. Cover the dentures with hand hot water (50°C). DO NOT USE BOILING WATER.
3. Add one tablet or sachet, close the lid, turn the switch to position II and let the machine vibrate ultra-sonically for 5 minutes.
4. Switch off after 5 minutes, rinse the dentures in cold water and brush lightly with a soft bristle brush. The same routine is used for all dentures, metal dentures and flexible dentures, mouth guards & gum shields

Place your mouth guard inside and use COLD WATER.

Turn to position II and let the machine vibrate ultra-sonically for 10 minutes. WARNING: Hot water can warp the gum shield.

Jewellery

- Place the jewellery on the tray. Cover the jewellery with hand hot water (50°C).
- Add five drops of washing up liquid and five drops of household ammonia solution and close the lid.
- Turn to position I and vibrate for 15 minutes. DO NOT USE PEARLS OR PEARL JEWELLERY AS YOU MAY DAMAGE THE PEARLS.
- The Swissedent Sonic Cleaner uses X3 AA batteries (NOT INCLUDED). To prevent the corrosion of the battery terminals in the damp atmosphere of the bathroom we recommend applying petroleum jelly to the terminals.

THIS LEAFLET IS PROVIDED FOR ADVICE ONLY.